

*Fast + Simple*  
**Focused on Veterinary Diagnostics**

# **FASTest® C. perfringens Toxin** ad us. vet.

## **Diarrhoea pathogen with high enterotoxic potential**

Fast test for the qualitative detection of *Clostridium perfringens* enterotoxin in feces of the dog, cat, goat and sheep lamb, calf, foal and piglet

### **Fast aetiological diagnostics**

### **Diarrhoea, haemorrhagical gastroenteritis**


### **Introduction of antibiogram (resistance avoidance)**

### **Initiation of specific therapy and prophylaxis**



- Simple test procedure with feces
- Fast test interpretation after 5 minutes
- Reliable clinical diagnostics
- Sensitivity 97.83 % & Specificity 98.15 %
- Storage at room temperature (15–25 °C)
- Long shelf life
- Compact test box with 2 or 10 tests

# ... **FASTest® C. perfringens Toxin** ad us. vet.

...  The gram positive anaerobic bacterium *Clostridium perfringens* belongs to the intestinal flora of many pets and farm animals and is facultative pathogenic. Inconvenient endogeneous (other basic diseases, diarrhoea pathogens, antibiotic therapies with massive reduction of intestinal flora etc.) and exogeneous (farming conditions, ex-treme changes of the food, stress etc.) factors can disturb the floral balance within the gut enabling *C. perfringens* to reproduce actively. Next to its ability to form extremely infectious and stable spores, the formation of lethal toxins is crucial for its pathogenicity.

The classification into types A–E is based on the different toxins that are produced. These toxins can cause extremely variable (mild to lethal progression forms) failures of the intestinal water and electrolyte balance in the different species like goat, sheep (e. g. dysentery of lambs: type B; pulpy kidney disease: type D), cattle (haemorrhagic enteritis: type A–E), foal (haemorrhagic necrotising enteritis: type A & C) and piglet (e. g. serous-catarhal enteritis: type A, necrotising enteritis: type C).

In the dog, especially serotype A occurs, producing 2 main toxins (toxin Alpha [α] and a Clostridia enterotoxin [CPE]), rarer serotype B (toxin Beta [β]). Both *C. perfringens* and its CPE can be detected also in healthy dog's feces. The CPE can be detected more often in dogs with diarrhoea compared to healthy dogs. For cats, to date reliable literature data concerning prevalence and clinical relevance are missing. Only by detection of *C. perfringens* in the feces, a disease caused by Clostridia is not diagnosable. Further investigation is necessary.

In a study in Switzerland, 54% of the *C. perfringens* isolates showed a reduced sensitivity towards metronidazole or 18% towards tetracycline. Because there is a general risk of resistance formation, it is recommended to identify the triggering pathogen in principle. By its high sensitivity and specificity, the use of **FASTest® C. perfringens Toxin** allows the veterinarian a rapid aetiological on-site diagnosis of a *C. perfringens* infection and subsequently the initiation of therapy as well as of necessary quarantine and prophylaxis measures.

## Test procedure



## Test interpretation



POSITIVE



NEGATIVE



Diarrhoea can have other reasons, too. Therefore, for diarrhoea symptoms it is generally advisable to use parallel tests **FASTest® BCV Strip**, **FASTest® C. diff 2T**, **FASTest® CCoV Strip**, **FASTest® CRYPTO Strip**, **FASTest® CRYPTO-ROTA D2T**, **FASTest® D4T bovine**, **FASTest® E.coli-K99 Strip**, **FASTest® GIARDIA Strip**, **FASTest® PARVO Card/Strip** or **FASTest® ROTA Strip** as well as **FASTest® CRP** canine or **FASTest® SAA** for the detection of an inflammation.

EN 01-2025



Distribution: