

## MultiCoder: Overview product carriers and accessories

Carriers and accessories		
<i>Single print</i>		
Name	Illustration	Ref.
Carrier for goblets: goblet Ø 35 mm		13145/0100
goblet Ø 65 mm		13145/0101
goblet hexagonal		13145/0102
Carrier for flat cassette		13145/0110
Carrier for petri dishes		13145/0300
Carrier for 5x cane marking labels		13145/0103



## Carriers and accessories

### Multiple print

Name	Illustration	Ref.
<p>Insert for multiple printing</p> <p>Basic requirement for multiple printing: combination of multiple print carrier and insert</p>		<p>13145/0200</p>
<p>Carrier for goblets:</p> <p>goblet Ø 7.1 mm (20x)</p> <p>goblet Ø 10 mm (15x)</p> <p>goblet Ø 13 mm (12x)</p>		<p>13145/0219</p> <p>13145/0220</p> <p>13145/0221</p>
<p>Carrier for straws:</p> <p>straws 0.25 ml (60x)</p> <p>straws 0.5 ml (36x)</p>		<p>13145/0201</p> <p>13145/0202</p>
<p>Adapter for IceCube racks:</p> <p>straws 0.5 ml (36x)</p>		<p>13145/0210</p>

## Carriers and accessories

### Accessories

Name	Illustration	Ref.
Loading block for 60 straws 0.25 ml		13145/0230
Loading block for 36 straws 0.5 ml		13145/0231

**Carrier – print object – allocation**

*Single print*

Product	Carrier/Ref.
Goblet, round, 35 mm 	 13145/0100
Goblet, round, 65 mm 	 13145/0101
Goblet, hexagonal 	 13145/0102
Flat cassette 	 13145/0110
Petri dishes 	 13145/0300

Cane marking labels



5 labels on one carrier  
print: 5 labels/cycle

13145/0103







## Carrier – print object – allocation










### Multiple print



Basic requirement for multiple printing: combination of multiple print carrier and insert

[13145/0200](#)

Product	Carrier/Ref.	Accessories/Ref.
Goblet, round, 7.1 mm 	 20 goblets on one carrier print: 1 goblet/cycle <a href="#">13145/0219</a>	
Goblet, round, 10 mm 	 15 goblets on one carrier print: 1 goblet/cycle <a href="#">13145/0220</a>	
Goblet, round, 13 mm 	 12 goblets on one carrier print: 1 goblet/cycle <a href="#">13145/0221</a>	

<p>Straw, 0.25 ml</p> 	 <p>60 straws on one carrier print: 4 straws/cycle</p> <p><a href="#">13145/0201</a></p>	 <p>Loading block straws 0.25 ml</p> <p><a href="#">13145/0230</a></p>
<p>Straw, 0.5 ml</p> 	 <p>36 straws on one carrier print: 3 straws/cycle</p> <p><a href="#">13145/0202</a></p>	 <p>Loading block straws 0.5 ml</p> <p><a href="#">13145/0231</a></p>
<p>Straws on IceCube racks</p> 	 <p>36 straws on one carrier print: 3 straws/cycle</p> <p><a href="#">13145/0210</a></p>	 <p>IceCube rack and loading block</p> <p><a href="#">16821/5036</a> &amp; <a href="#">16821/5037</a></p>